



ENFit™  
Connectors


# Your Clear Connection to Enteral Safety: Reducing the Risk of Misconnections

Cross Reference Product Catalog for ENFit™ Connector Transition Assistance



*Issued: July 2016*





## Key to Product Code Transition to ENFit™

Product Category	Old Product Code Scheme	New Product Code Scheme
NG/NI Tubes, AIV2	20-xxxx, 20-XXXXAIV2	40-xxxx, 42-xxxx (Sterile)
CORTRAK 2 EAS Feeding Tubes	20-xxxxTRAK2	40-xxxxTRAK2
CORTRAK/CORGRIP Convenience Kits	25-xxxxTRAK2G	45-xxxxTRAK2G
NICU AIV → CORLITE	20-xxxxAIV	41-xxxx
ENDO/F → CORSCOPE For Endoscopic Placement	30-xxxx	50-xxxx
FARRELL Valve System	20-4100	43-4100
PEG Tubes and Accessories	30-xxxx	50-xxxx
Jejunostomy (J-Tube) Extension Tubes (for use with MAX PEG Tubes)	30-xxxx	50-xxxx
Gastrostomy Tubes (G-Tubes)		
	Dual 32-xxxx	52-xxxx
	Triple 33-xxxx	53-xxxx
cuBBy Low Profile G-Tubes (LPGD)	35-xxxx	55-xxxx
LPGD Accessories	35-xxxx	56-xxxx

## ENFit™ Transition Connectors—No ENFit™ Codes Needed

ENFIT™ CATALOG #	PRODUCT	PACKAGED
F00071	<b>Step (“Christmas Tree”) Connector:</b> Allows an ENFit feeding administration set or enteral gastric decompression device to connect to a standard/non-ENFit feeding tube	30 per case
F00073	<b>Oral Syringe Connector:</b> Allows an Oral Syringe to be connected to an ENFit feeding tube	30 per case
F00074	<b>Transition Funnel Connector:</b> Allows an ENFit feeding tube to accept a Step (“Christmas Tree”) Connector or a Catheter Tip Syringe	30 per case

# CORTRAK<sup>®</sup> 2 Enteral Access System

## CORTRAK<sup>®</sup> 2 Enteral Access System (EAS<sup>™</sup>)—No ENFit<sup>™</sup> Codes Needed

CATALOG #	PRODUCT	PACKAGED
20-0950	CORTRAK <sup>®</sup> 2 Enteral Access System (EAS <sup>™</sup> ) All-In-One Monitor, Integrated Visual Display Terminal (VDT), Touch Screen and Embedded Computing System	1 per case
	Smart Receiver Unit	Include
	Power Supply	Included
	Power Cord	Included
	Interconnect Cable	Included
	Operator's Manual with Service Guide	Included
	Tip Guide Booklet	Included
20-0951	CORTRAK 2 EAS Accessory Kit	1 per case

## CORTRAK<sup>®</sup> 2 Enteral Access System (EAS<sup>™</sup>) Feeding Tubes

CORTRAK<sup>®</sup> 2 Enteral Access System Electromagnetic Transmitting Stylet and CORFLO<sup>®</sup> Feeding Tube; ULTRA, non-weighted, anti-clog bolus.

OLD CATALOG #	ENFIT <sup>™</sup> CATALOG #	PRODUCT	PACKAGED
20-9368TRAK2	40-9368TRAK2	8 FR 36" Length	10 tubes per case
20-9361TRAK2	40-9361TRAK2	10 FR 36" Length	10 tubes per case
20-9438TRAK2	40-9438TRAK2	8 FR 43" Length	10 tubes per case
20-9432TRAK2	40-9432TRAK2	12 FR 43" Length	10 tubes per case
20-9431TRAK2	40-9431TRAK2	10 FR 43" Length	10 tubes per case
20-9558TRAK2	40-9558TRAK2	8 FR 55" Length	10 tubes per case
20-9551TRAK2	40-9551TRAK2	10 FR 55" Length	10 tubes per case

## CORTRAK<sup>®</sup> 2 Enteral Access System (EAS<sup>™</sup>) Feeding Tubes with Anti-IV<sup>®</sup> (AIV) Connector

*Anti-IV<sup>®</sup> CORTRAK Feeding Tubes are being phased out as ENFit<sup>™</sup> equivalent products are introduced (CORLITE<sup>™</sup>).*

CORTRAK<sup>®</sup> 2 EAS<sup>™</sup> Electromagnetic Transmitting Stylet and CORFLO<sup>®</sup> Feeding Tubes; ULTRA-Lite, Non-weighted, Anti-clog bolus.

OLD CATALOG #	ENFIT <sup>™</sup> CATALOG #	PRODUCT	PACKAGED
20-9368TRAK2A	40-9368TRAK2	8 FR, 36" Length	10 tubes per case
20-9361TRAK2A	40-9361TRAK2	10 FR, 36" Length	10 tubes per case
20-9438TRAK2A	40-9438TRAK2	8 FR, 43" Length	10 tubes per case
20-9432TRAK2A	40-9432TRAK2	12 FR, 43" Length	10 tubes per case
20-9431TRAK2A	40-9431TRAK2	10 FR, 43" Length	10 tubes per case
20-9558TRAK2A	40-9558TRAK2	8 FR, 55" Length	10 tubes per case
20-9551TRAK2A	40-9551TRAK2	10 FR, 55" Length	10 tubes per case

**CORTRAK® 2 Enteral Access System (EAS™) Replacement Parts—  
No ENFit™ Codes Needed**

CATALOG #	PRODUCT	PACKAGED
20-0952	Smart Receiver Unit	1 per case
20-0966	Interconnect Cable	1 per case
20-0953	Li-Ion Battery Pack	1 per case

**CORTRAK® 2 Enteral Access System (EAS™) Accessory Items—  
No ENFit™ Codes Needed**

CATALOG #	PRODUCT	PACKAGED
20-0918	QL220 Plus Printer (includes 2 rechargeable batteries)	1 per case
20-0919	QL220 Plus Quad Charger	1 per case
20-0911	Printer Paper Rolls (125 per Roll)	3 per case
20-0916	Printer Replacement Battery (QL220 Plus)	1 per case
20-0920	Stand (Includes SRU/RU Cover Holder)	1 per case
20-0930	Plastic Covers (Sleeves) for SRU/RU (4 rolls, 25 bags per roll)	100 per case
20-0955	Stand Bracket Kit	1 per case



**NAVIGATØR® BioNavigation™ System**

CATALOG #	PRODUCT	PACKAGED
60-3020	NAVIGATØR® Instrument with Case and Extension Cable	1 per case
60-2018	1.1 FR MAPCath® Sensor Stylet, 36 mm, 106 cm	10 per case
60-2011	1.4 FR MAPCath® Sensor Stylet, 46 mm, 106 cm	10 per case
60-3002	NAVIGATØR® Extension Cable	1 per case
60-2048	NAVIGATØR® Sterility Sleeve, 6" x 48"	25 per case
60-2096	NAVIGATØR® Sterility Sleeve, 6" x 96"	25 per case



## NASOENTERIC

Ask for the **YELLOW Tube...**

Opaque polyurethane feeding tube (unless otherwise noted), medication/flushing port, pre-inserted stylet, tungsten weighted tip, C-19™ hydrophilic lubricant on tube tip and internal lumen, radiopaque

### CORFLO®-ULTRA Pill NG—With Stylet

Anti-clog, pill shape bolus, 3 gm\* adult weighted, with stylet

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-6438	40-6438	8 FR 43" Length	10 tubes per case
20-6431	40-6431	10 FR 43" Length	10 tubes per case



\*Approximate weight

### CORFLO®-ULTRA NG—With Stylet

Anti-clog, smooth shape bolus, 3 gm\* adult weighted, with stylet

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-7366	40-7366	6 FR 36" Length, 1.5 gm*	10 tubes per case
20-7368	40-7368	8 FR 36" Length	10 tubes per case
20-7361	40-7361	10 FR 36" Length	10 tubes per case
20-7438	40-7438	8 FR 43" Length	10 tubes per case
20-7431	40-7431	10 FR 43" Length	10 tubes per case
20-7432	40-7432	12 FR 43" Length	10 tubes per case



\*Approximate weight

### CORFLO®-ULTRA NG

Anti-clog, smooth shape bolus, 3 gm\* adult weighted, without stylet

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-5361	40-5361	10 FR 36" Length	10 tubes per case
20-5362	40-5362	12 FR 36" Length	10 tubes per case
20-5438	40-5438	8 FR 43" Length	10 tubes per case
20-5431	40-5431	10 FR 43" Length	10 tubes per case
20-5432	40-5432	12 FR 43" Length	10 tubes per case



\*Approximate weight

### CORFLO®-ULTRA 7 NG—With Stylet

Anti-clog, smooth shape bolus, 7 gm\* adult weighted, with stylet

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-4438	40-4438	8 FR 43" Length	10 tubes per case
20-4431	40-4431	10 FR 43" Length	10 tubes per case
20-4432	40-4432	12 FR 43" Length	10 tubes per case



\*Approximate weight

## CORFLO®-ULTRA Ped NG—With Stylet

Anti-clog, smooth shape bolus, pediatric weighted (1gm\*), with stylet

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-8226	40-8226	6 FR 22" Length	10 tubes per case
20-8366	40-8366	6 FR 36" Length	10 tubes per case
20-8368	40-8368	8 FR 36" Length	10 tubes per case



\*Approximate weight

## CORFLO®-ULTRA Lite NG—With Stylet

Anti-clog, bowed tip bolus, non-weighted, adult-pediatric-neonatal, with stylet

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-9225*	42-9225*	5 FR 22" Length	10 tubes per case
20-9226	40-9226	6 FR 22" Length	10 tubes per case
20-9228	40-9228	8 FR 22" Length	10 tubes per case
20-9366	40-9366	6 FR 36" Length	10 tubes per case
20-9368	40-9368	8 FR 36" Length	10 tubes per case
20-9361	40-9361	10 FR 36" Length	10 tubes per case
20-9438	40-9438	8 FR 43" Length	10 tubes per case
20-9431	40-9431	10 FR 43" Length	10 tubes per case
20-9432	40-9432	12 FR 43" Length	10 tubes per case



\*STERILE

## CORFLO®-ULTRA Lite NG

Anti-clog, bowed tip bolus, non-weighted, adult-pediatric-neonatal, without stylet

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-1155*	42-1155*	5 FR 15" Length	10 tubes per case
20-1156*	42-1156*	6 FR 15" Length	10 tubes per case
N/A	42-1158*	8 FR 15" Length	10 tubes per case
N/A	42-1225*	5 FR 22" Length	10 tubes per case
20-1226	42-1226*	6 FR 22" Length	10 tubes per case
N/A	42-1228*	8 FR 22" Length	10 tubes per case
20-1365*	42-1365*	5 FR 36" Length	10 tubes per case
20-1366	40-1366	6 FR 36" Length	10 tubes per case
20-1368	40-1368	8 FR 36" Length	10 tubes per case
20-1361	40-1361	10 FR 36" Length	10 tubes per case
20-1362	40-1362	12 FR 36" Length	10 tubes per case
20-1438	40-1438	8 FR 43" Length	10 tubes per case
20-1431	40-1431	10 FR 43" Length	10 tubes per case
20-1432	40-1432	12 FR 43" Length	10 tubes per case



\*STERILE

## CORFLO®-ULTRA Lite Clear NG

This category is being replaced by either CORLITE™ (Single access port, clear with radiopaque stripe: 41-xxxx) or CORFLO® ULTRA Lite tubes (Y access port, opaque Yellow tubes, ULTRA Anti-Clog Feeding Port: 42-1xxx).

OLD CATALOG #	CORLITE™ ENFIT™ CATALOG #	CORFLO® ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-C1225*	41-1225*	42-1225*	5 FR 22" Length	10 tubes per case
20-C1226*	41-1226*	42-1226*	6 FR 22" Length	10 tubes per case
20-C1228*	41-1228*	42-1228*	8 FR 22" Length	10 tubes per case



\*STERILE

## CORFLO® – Controller Enteral Feeding Tubes

### CORFLO®-Controller Pill NG—With Stylet

Pill shape bolus, 3 gm\* adult weighted, with stylet

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-2438	40-2438	8 FR 43" Length	10 tubes per case
20-2431	40-2431	10 FR 43" Length	10 tubes per case
20-2551	40-2551	10 FR 55" Length	10 tubes per case

\*Approximate weight

### CORFLO®-Controller NG—With Stylet

Smooth shape bolus, 3 gm\* adult weighted, with stylet

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-3438	40-3438	8 FR 43" Length	10 tubes per case
20-3431	40-3431	10 FR 43" Length	10 tubes per case
20-3551	40-3551	10 FR 55" Length	10 tubes per case

\*Approximate weight

### CORFLO®-Controller 7 NG—With Stylet

7 gm\* adult weighted, with stylet

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-0438	40-0438	8 FR 43" Length	10 tubes per case
20-0431	40-0431	10 FR 43" Length	10 tubes per case
20-0432	40-0432	12 FR 43" Length	10 tubes per case

\*Approximate weight

# CORFLO®-CORSCOPE® Nasointestinal Endoscopically Placed Feeding Tube

(Previously called ENDO/F)

60" CORFLO®-ULTRA Anti-Clog Feeding Tube, 300 cm Teflon Coated Flexible Guide Wire (0.035" diameter) in dispenser sheath, and 10 FR, 22" length Nasal/Oral Transfer Tube

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
30-4602	50-4602	CORFLO®-ENDO/F NI 12 FR 7 gm* weighted	1 per case
30-9601	50-9601	CORFLO®-ENDO/F NI 10 FR Non-weighted	1 per case



\*Approximate weight

## Specialty Small Bowel Feeding (SBF) Tubes

SBF Convenience Kits are being eliminated with the introduction of ENFit™. The same tube with ENFit™ is still available as a “non-K” (Kit) code as per below.

For those interested in bedside placement of small bowel feeding tubes, please contact your sales representative to learn about the CORTRAK® Enteral Access System.

Anti-clog, bowed tip bolus, non-weighted, with stylet

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-9431K	40-9431	10 FR 43" Length	10 kits per case
20-9551K	40-9551	10 FR 55" Length	10 kits per case

### CORFLO®-ULTRA Lite SBF NG—With Stylet

Anti-clog, bowed tip bolus, non-weighted adult tube, with stylet

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-9558	40-9558	8 FR 55" Length	10 tubes per case
20-9551	40-9551	10 FR 55" Length	10 tubes per case



### CORFLO®-Controller NG—With Stylet

Smooth shape bolus, 3 gm\* adult weighted, with stylet

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-3551	40-3551	10 FR 55" Length	10 tubes per case

\*Approximate weight

### CORFLO®-Controller Pill NG—With Stylet

Pill shape bolus, 3 gm\* adult weighted, with stylet

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-2551	40-2551	10 FR 55" Length	10 tubes per case

\*Approximate weight





## CORFLO® Anti-IV® (AIV2) Feeding Tubes

CORFLO® Anti-IV® (AIV2) Feeding Tubes are being phased out as ENFit™ equivalent products are introduced.

The following feeding tubes are designed to be used with standard enteral pump administration sets.

### CORFLO ULTRA-Lite NG—With Stylet and Anti-IV® (AIV) Connectors

Anti-clog bolus, non-weighted, adult-pediatric-neonatal with stylet

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-9225AIV2*	42-9225*	5 FR, 22" Length	10 tubes per case
20-9226AIV2	40-9226	6 FR, 22" Length	10 tubes per case
20-9228AIV2	40-9228	8 FR, 22" Length	10 tubes per case
20-9366AIV2	40-9366	6 FR, 36" Length	10 tubes per case
20-9368AIV2	40-9368	8 FR, 36" Length	10 tubes per case
20-9361AIV2	40-9361	10 FR, 36" Length	10 tubes per case
20-9438AIV2	40-9438	8 FR, 43" Length	10 tubes per case
20-9431AIV2	40-9431	10FR, 43" Length	10 tubes per case
20-9432AIV2	40-9432	12 FR, 43" Length	10 tubes per case
20-9558AIV2	40-9558	8 FR, 55" Length	10 tubes per case
20-9551AIV2	40-9551	10FR, 55" Length	10 tubes per case



\*STERILE

### CORFLO® ULTRA-Lite NG Without Stylet and With Anti-IV® (AIV) Connectors

Anti-clog bolus, non-weighted, adult-pediatric-neonatal without stylet

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-1155AIV2*	42-1155*	5 FR, 15" Length	10 tubes per case
20-1156AIV2*	42-1156*	6 FR, 15" Length	10 tubes per case
20-1226AIV2	42-1226*	6 FR, 22" Length	10 tubes per case
20-1365AIV2*	42-1365*	5 FR, 36" Length	10 tubes per case
20-1366AIV2	40-1366	6 FR, 36" Length	10 tubes per case
20-1368AIV2	40-1368	8 FR, 36" Length	10 tubes per case
20-1361AIV2	40-1361	10 FR, 36" Length	10 tubes per case
20-1362AIV2	40-1362	12 FR, 36" Length	10 tubes per case
20-1438AIV2	40-1438	8 FR, 43" Length	10 tubes per case
20-1431AIV2	40-1431	10FR, 43" Length	10 tubes per case
20-1432AIV2	40-1432	12 FR, 43" Length	10 tubes per case



\*STERILE

## CORFLO®-ULTRA NG with Stylet and Anti-IV® (AIV) Connectors

Anti-clog bolus, adult weight (3gm\*) with stylet

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-7366AIV2	40-7366	6 FR, 36" Length	10 tubes per case
20-7368AIV2	40-7368	8 FR, 36" Length	10 tubes per case
20-7361AIV2	40-7361	10 FR, 36" Length	10 tubes per case
20-7438AIV2	40-7438	8 FR, 43" Length	10 tubes per case
20-7431AIV2	40-7431	10 FR, 43" Length	10 tubes per case
20-7432AIV2	40-7432	12 FR, 43" Length	10 tubes per case



\*Approximate weight

## CORFLO®-ULTRA PED NG with Stylet and Anti-IV® (AIV) Connectors

Anti-clog bolus, pediatric weight (1gm\*) with stylet

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-8226AIV2	40-8226	6 FR, 22" Length	10 tubes per case
20-8366AIV2	40-8366	6 FR, 36" Length	10 tubes per case
20-8368AIV2	40-8368	8 FR, 36" Length	10 tubes per case



\*Approximate weight



## CORFLO® Neonatal Anti-IV® (AIV)

*CORPAK's Neonatal AIV feeding tubes will now be available as CORLITE™*

Single ENFit™ Feeding Port, Clear polyurethane feeding tube with radiopaque stripe, Sterile

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-1155AIV	41-1155	5 FR, 15" Length	10 tubes per case
20-1156AIV	41-1156	6 FR, 15" Length	10 tubes per case
20-1158AIV	41-1158	8 FR, 15" Length	10 tubes per case
20-1225AIV	41-1225	5 FR, 22" Length	10 tubes per case
20-1226AIV	41-1226	6 FR, 22" Length	10 tubes per case
20-1228AIV	41-1228	8 FR, 22" Length	10 tubes per case
20-1365AIV	41-1365	5 FR, 36" Length	10 tubes per case
20-1366AIV	41-1366	6 FR, 36" Length	10 tubes per case
20-1368AIV	41-1368	8 FR, 36" Length	10 tubes per case

\*Approximate weight

## Neonatal Extension Sets

Clear DEHP Free microbore tubing, ENFit™ Connectors, Sterile

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-1030AIV	41-1030	30" Extension Set	50 per case
N/A	41-1030Y	30" Extension Set with in-line Y/Med Port	50 per case
20-1060AIV	41-1060	60" Extension Set	50 per case
N/A	41-1060Y	60" Extension Set with in-line Y/Med Port	50 per case
20-1060AIV-S	41-1060 or 41-1060Y	60" Extension Set or Extension Set with in-line Y/Med Port	50 per case

## Enteral Accessory Products



### FARRELL® Valve System with Equilibrium Technology™

Leading gastric pressure device, designed to help those who suffer from poor gastric motility, pain, and bloating. Maximizes enteral feeding by providing a channel to constantly decompress the stomach, allowing the stomach to feed at its own pace.

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-4100	43-4100	FARRELL® Valve Enteral Gastric Pressure Relief System	30 per case



### CORGRIP™ Tube Retention System—No ENFit™ Codes Needed

CATALOG #	PRODUCT	PACKAGED
25-008	For use with 8 FR NG/NI feeding tubes	10 per case
25-010	For use with 10 FR NG/NI feeding tubes	10 per case
25-012	For use with 12 FR NG/NI feeding tubes	10 per case

### CORTRAK®/CORGRIP™ Tube Retention System Convenience Kits

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
25-9368TRAK2G	45-9368TRAK2G	36" 8 FR CORTRAK®/CORGRIP™ Convenience Kit	5 per case
25-9361TRAK2G	45-9361TRAK2G	36" 10 FR CORTRAK®/CORGRIP™ Convenience Kit	5 per case
25-9438TRAK2G	45-9438TRAK2G	43" 8 FR CORTRAK®/CORGRIP™ Convenience Kit	5 per case
25-9431TRAK2G	45-9431TRAK2G	43" 10 FR CORTRAK®/CORGRIP™ Convenience Kit	5 per case
25-9432TRAK2G	45-9432TRAK2G	43" 12 FR CORTRAK®/CORGRIP™ Convenience Kit	5 per case
25-9558TRAK2G	45-9558TRAK2G	55" 8 FR CORTRAK®/CORGRIP™ Convenience Kit	5 per case
25-9551TRAK2G	45-9551TRAK2G	55" 10 FR CORTRAK®/CORGRIP™ Convenience Kit	5 per case



In-Patient Tube Clearing System

### TubeClear® Control Box—No ENFit™ Codes Needed

CATALOG NUMBER	PRODUCT	PACKAGED
101-US	Control Box—US	1 each

### TubeClear® Clearing Stems for NG/NI Feeding or Sump Tubes—No ENFit™ Codes Needed

CATALOG #	PRODUCT	COLOR CODE	PACKAGED
5-NE-1036	36" Length; 10–18 FR	Yellow	5 per case
5-NE-1042	42" Length; 10–18 FR	Gray	5 per case
5-NE-1043	43" Length; 10–18 FR	Purple	5 per case
5-NE-1045	45" Length; 10–18 FR	Orange	5 per case
5-NE-1048	48" Length; 10–18 FR	Clear	5 per case
5-NE-1050	50" Length; 10–18 FR	Blue	5 per case
5-NE-1055	55" Length; 10–18 FR	Green	5 per case

### TubeClear® Clearing Stems for PEG, Gastrostomy, or Jejunostomy Tubes—No ENFit™ Codes Needed

CATALOG #	PRODUCT	COLOR CODE	PACKAGED
5-G-1008	8" Length; 10–18 FR	Yellow	5 per case
5-G-1009	9" Length; 10–18 FR	Gray	5 per case
5-G-1010	10" Length; 10–18 FR	Purple	5 per case
5-G-1011	11" Length; 10–18 FR	Orange	5 per case
5-G-1012	12" Length; 10–18 FR	Clear	5 per case
5-G-1014	14" Length; 10–18 FR	Blue	5 per case

### Other Enteral Accessory Products

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
20-0002	No Change	CLOG ZAPPER® Enteral Feeding Tube Declogging System	10 per case
20-1007	40-1007	Enteral Y Extension Set	30 per case

## CORFLO®-MAX PEG Kits & Accessory Products

### Pull Technique

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	FR SIZE/BOLSTER DESIGN	PACKAGED
30-5020	50-5020	CORFLO®-MAX PEG Kit (Contains Lidocaine & Povidone Iodine)	20 FR—Conical	2 per case
30-5020-S5	50-5020-S5	CORFLO®-MAX PEG Kit with 5 FR Snare	20 FR—Conical	2 per case
30-5020-S7	50-5020-S7	CORFLO®-MAX PEG Kit (Contains Lidocaine & Povidone Iodine)	20 FR—Conical	2 per case
30-5320	50-5320	CORFLO®-MAX SAFETY PEG Kit	20 FR—Conical	2 per case
30-5320-S7	50-5320-S7	CORFLO®-MAX SAFETY PEG Kit with 7 FR Snare	20 FR—Conical	2 per case
30-5520	50-5520	CORFLO®-MAX PEG Kit—Basic Components Tray (Contains No Medicines)	20 FR—Conical	2 per case
30-6020	50-6020	CORFLO®-MAX PEG Kit (Contains Lidocaine & Povidone Iodine)	20 FR—Ring	2 per case
30-6020-S7	50-6020-S7	CORFLO®-MAX PEG Kit with 7 FR Snare	20 FR—Ring	2 per case
30-6016	50-6016	CORFLO®-MAX PEG Kit (Contains Lidocaine & Povidone Iodine)	16 FR—Ring	2 per case
30-6016-S5	50-6016-S5	CORFLO®-MAX PEG Kit with 5 FR Snare (Contains Lidocaine & Povidone Iodine)	16 FR—Ring	2 per case
30-6012	50-6012	CORFLO®-MAX PEG Kit (Contains Lidocaine & Povidone Iodine)	12 FR—Ring	2 per case
30-6012-S5	50-6012-S5	CORFLO®-MAX PEG Kit (Contains Lidocaine & Povidone Iodine)	12 FR—Ring	2 per case
30-6316	50-6316	CORFLO®-MAX SAFETY PEG Kit	16 FR—Ring	2 per case
30-6316-S5	50-6316-S5	CORFLO®-MAX SAFETY PEG Kit with 5 FR Snare	16 FR—Ring	2 per case
30-6312	50-6312	CORFLO®-MAX SAFETY PEG Kit	12 FR—Ring	2 per case
30-6312-S5	50-6312-S5	CORFLO®-MAX SAFETY PEG Kit with 5 FR Snare	12 FR—Ring	2 per case
30-6516	50-6516	CORFLO®-MAX PEG Kit— Basic Components Tray	16 FR—Ring	2 per case
30-6512	50-6512	CORFLO®-MAX PEG Kit— Basic Components Tray	12 FR—Ring	2 per case
30-6520	50-6520	CORFLO®-MAX PEG Kit— Basic Components Tray	20 FR—Ring	2 per case

## Push Technique

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	FR SIZE/BOLSTER DESIGN	PACKAGED
30-4020	50-4020	CORFLO®-MAX PEG Kit (Contains Lidocaine & Povidone Iodine)	20 FR—Conical	2 per case
30-4020-S5	50-4020-S5	CORFLO®-MAX PEG Kit with 5 FR Snare (Contains Lidocaine & Povidone Iodine)	20 FR—Conical	2 per case
30-4020-S7	50-4020-S7	CORFLO®-MAX PEG Kit with 7 FR Snare (Contains Lidocaine & Povidone Iodine)	20 FR—Conical	2 per case
30-4320	50-4320	CORFLO® Safety MAX PEG Kit	20 FR—Conical	2 per case
30-4320-S7	50-4320-S7	CORFLO® Safety MAX PEG Kit with 7 FR Snare	20 FR—Conical	2 per case
30-7020	50-7020	CORFLO®-MAX PEG Kit (Contains Lidocaine & Povidone Iodine)	20 FR—Ring	2 per case
30-7020-S7	50-7020-S7	CORFLO®-MAX PEG Kit with 7 FR Snare	20 FR—Ring	2 per case

## CORLOCK®-CORPORT® Y-Adapter Repair Kits

### For Use With CORFLO®-MAX PEG Tubes

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
30-6120	50-6120	CORLOCK®-CORPORT® Y-Adapter Repair Kit—20 FR	5 per case
30-6116	50-6116	CORLOCK®-CORPORT® Y-Adapter Repair Kit—16 FR	5 per case
30-6112	50-6112	CORLOCK®-CORPORT® Y-Adapter Repair Kit—12 FR	5 per case
30-6220	50-6220	CORLOCK®-TWOOMEY Y-Adapter Repair Kit—20 FR	5 per case
30-6216	50-6216	CORLOCK®-TWOOMEY Y-Adapter Repair Kit—16 FR	5 per case

## CORFLO®-ULTRA Jejunal Tubes

### For Use With CORFLO®-MAX PEG Tubes

Weighted Jejunal Tube with CORFLO®-ULTRA Anti-Clog Feeding Port, Weighted, with Grip-Tip, Adjust-a-Sleeve, CORLOCK®-TWOOMEY Y-Adapter, Silicone Plug, and Stylet with Cap

OLD CATALOG #	ENFIT™ CATALOG #	PRODUCT	PACKAGED
30-7366	50-7366	6 FR, 36"; 3 gm* Weighted For Use With 16 FR PEG Tube	5 per case
30-7368	50-7368	8 FR, 36"; 3 gm* Weighted For Use With 20 FR PEG Tube	5 per case
30-7361	50-7361	10 FR, 36"; 3 gm* Weighted For Use With 20 FR PEG Tube	5 per case



\*Approximate weight

# CORFLO<sup>®</sup>-ULTRA Over-the-Wire Jejunal Tubes

## For Use With 20 FR CORFLO<sup>®</sup>-MAX PEG Tubes

Non-Weighted Jejunal Tube with CORFLO<sup>®</sup>-ULTRA Anti-Clog Feeding Port and Adjust-a-Sleeve, CORLOCK<sup>®</sup> - TWOOOMEY Y-Adapter, 260 cm Teflon Coated Flexible Guide Wire (0.035" diameter) in Dispenser Sheath and Silicone Plug.

OLD CATALOG #	ENFIT <sup>™</sup> CATALOG #	PRODUCT	PACKAGED
30-1438	50-1438	8 FR, 43" Length	2 per case
30-1431	50-1431	10 FR, 43" Length	2 per case



## CORFLO<sup>®</sup>-Dual Gastrostomy Tube Kits GT

Clear silicone gastrostomy tubes, inflatable balloon, retention bolster, safety plug, 9" length, packaged with prefilled syringe, lube jelly, and gauze. STERILE

OLD CATALOG #	ENFIT <sup>™</sup> CATALOG #	PRODUCT	PACKAGED
32-0512	52-0512	12 FR, 5 cc Balloon	2 per case
32-0514	52-0514	14 FR, 5 cc Balloon	2 per case
32-0516	52-0516	16 FR, 5 cc Balloon	2 per case
32-0518	52-0518	18 FR, 5 cc Balloon	2 per case
32-2016	52-2016	16 FR, 20 cc Balloon	2 per case
32-2018	52-2018	18 FR, 20 cc Balloon	2 per case
32-2020	52-2020	20 FR, 20 cc Balloon	2 per case
32-2022	52-2022	22 FR, 20 cc Balloon	2 per case
32-2024	52-2024	24 FR, 20 cc Balloon	2 per case

## CORFLO<sup>®</sup>-Triple Gastrostomy Tube Kits GT

Clear silicone gastrostomy tubes, inflatable balloon, retention bolster, safety plug, 9" length, packaged with prefilled syringe, lube jelly, gauze. STERILE

OLD CATALOG #	ENFIT <sup>™</sup> CATALOG #	PRODUCT	PACKAGED
33-2016	53-2016	16 FR, 20 cc Balloon	2 per case
33-2018	53-2018	18 FR, 20 cc Balloon	2 per case
33-2020	53-2020	20 FR, 20 cc Balloon	2 per case
33-2022	53-2022	22 FR, 20 cc Balloon	2 per case
33-2024	53-2024	24 FR, 20 cc Balloon	2 per case

# CORFLO-cuBBy®

## Low Profile Gastrostomy Device (LPGD) Kit

Clear silicone Low Profile Gastrostomy Device in tray-STERILE, with non-sterile Care Giver Kit

OLD CATALOG #	†ENFIT™ CATALOG #	PRODUCT	PACKAGED
35-1210	55-1210	12 FR, x 1.0 cm	1 per case
35-1215	55-1215	12 FR, x 1.5 cm	1 per case
35-1220	55-1220	12 FR, x 2.0 cm	1 per case
35-1225	55-1225	12 FR, x 2.5 cm	1 per case
35-1230	55-1230	12 FR, x 3.0 cm	1 per case
35-1235	55-1235	12 FR, x 3.5 cm	1 per case
35-1410	55-1410	14 FR, x 1.0 cm	1 per case
35-1415	55-1415	14 FR, x 1.5 cm	1 per case
35-1420	55-1420	14 FR, x 2.0 cm	1 per case
35-1425	55-1425	14 FR, x 2.5 cm	1 per case
35-1430	55-1430	14 FR, x 3.0 cm	1 per case
35-1435	55-1435	14 FR, x 3.5 cm	1 per case
35-1440	55-1440	14 FR, x 4.0 cm	1 per case
35-1445	55-1445	14 FR, x 4.5 cm	1 per case
35-1610	55-1610	16 FR, x 1.0 cm	1 per case
35-1615	55-1615	16 FR, x 1.5 cm	1 per case
35-1620	55-1620	16 FR, x 2.0 cm	1 per case
35-1625	55-1625	16 FR, x 2.5 cm	1 per case
35-1630	55-1630	16 FR, x 3.0 cm	1 per case
35-1635	55-1635	16 FR, x 3.5 cm	1 per case
35-1640	55-1640	16 FR, x 4.0 cm	1 per case
35-1645	55-1645	16 FR, x 4.5 cm	1 per case
35-1810	55-1810	18 FR, x 1.0 cm	1 per case

OLD CATALOG #	†ENFIT™ CATALOG #	PRODUCT	PACKAGED
35-1815	55-1815	18 FR, x 1.5 cm	1 per case
35-1820	55-1820	18 FR, x 2.0 cm	1 per case
35-1825	55-1825	18 FR, x 2.5 cm	1 per case
35-1830	55-1830	18 FR, x 3.0 cm	1 per case
35-1835	55-1835	18 FR, x 3.5 cm	1 per case
35-1840	55-1840	18 FR, x 4.0 cm	1 per case
35-1845	55-1845	18 FR, x 4.5 cm	1 per case
35-2010	55-2010	20 FR, x 1.0 cm	1 per case
35-2015	55-2015	20 FR, x 1.5 cm	1 per case
35-2020	55-2020	20 FR, x 2.0 cm	1 per case
35-2025	55-2025	20 FR, x 2.5 cm	1 per case
35-2030	55-2030	20 FR, x 3.0 cm	1 per case
35-2035	55-2035	20 FR, x 3.5 cm	1 per case
35-2040	55-2040	20 FR, x 4.0 cm	1 per case
35-2045	55-2045	20 FR, x 4.5 cm	1 per case
35-2410	55-2410	24 FR, x 1.0 cm	1 per case
35-2415	55-2415	24 FR, x 1.5 cm	1 per case
35-2420	55-2420	24 FR, x 2.0 cm	1 per case
35-2425	55-2425	24 FR, x 2.5 cm	1 per case
35-2430	55-2430	24 FR, x 3.0 cm	1 per case
35-2435	55-2435	24 FR, x 3.5 cm	1 per case
35-2440	55-2440	24 FR, x 4.0 cm	1 per case
35-2445	55-2445	24 FR, x 4.5 cm	1 per case

## CORFLO-cuBBy® Accessories

OLD CATALOG #	†ENFIT™ CATALOG #	PRODUCT	PACKAGED
35-9024	56-9024	Right Angle Feeding Set for use with the CORFLO-cuBBy® (24" length)	5 per case
35-9012	56-9012	Right Angle Feeding Set for use with the CORFLO-cuBBy® (12" length)	5 per case
35-0112	56-0112	Bolus Feeding Set for use with the CORFLO-cuBBy®	5 per case
31-5555	no change	Stoma Measuring Device for use with the CORFLO-cuBBy®	5 per case

NOTE: Product availability subject to change without notice.

†Pending FDA Clearance





## Lifenest<sup>®</sup> Sleep System

### Lifenest<sup>®</sup> Neonatal Mattress

Lifenest Sleep System is a low pressure mattress with a unique woven cradel offering diverse applications for the Neonatal (NICU) setting.

CATALOG #	PRODUCT	DIMENSIONS	COMPATIBILITY	PACKAGED
HMINC001UBIA	Incubator Mattress	63 x 48 x 6 cm; 24.8" x 18.9" x 2.3"	GE Giraffe™; Dräger Caleo™	6 per box
HMINC002UBIA	Incubator Mattress	75 x 39.5 x 6 cm; 29.5" x 15.5" x 2.3"	Dräger Isolette™ 8000	6 per box
HMINC003UBIA	Incubator Mattress	74 x 36 x 6 cm; 29.1" x 14.2" x 2.3"	ATOM V2100-G; ATOM V2200	6 per box
HMINC004UBIA	Incubator Mattress	61 x 40 x 6 cm; 24" x 15.7" x 2.3"	NATUS NatalCare™ LX	6 per box
HMINC005UBIA	Incubator Mattress	63.5 x 35 x 6 cm; 25" x 13.8" x 2.3"	NATUS NatalCare™ PC-305	6 per box
HMBAS001UBIA	Bassinet Mattress	66 x 31.5 x 6 cm; 26" x 12.4" x 2.3"	Most Hospital Bassinets	6 per box
HSPM004UBIA	Crib Mattress	85 x 68 x 10 cm; 33.5" x 26.8" x 4"	Most Hospital Cribs	2 per box
HSPM005UBIA	Crib Mattress	85 x 59 x 10 cm; 33.5" x 23.2" x 4"	Most Hospital Cribs	2 per box

### Lifenest<sup>®</sup> Disposable Sheets

Lifenest breathable sheets enhance the performance of the Lifenest mattress and are available in multiple sizes to accommodate any hospital mattress.

CATALOG #	PRODUCT	DIMENSIONS	COMPATIBILITY	PACKAGED
HC1MINC001UBIA	Incubator Sheet	63 x 48 x 6 cm; 24.8" x 18.9" x 2.3"	GE Giraffe™; Dräger Caleo™	100 per box
HC2MINC001UBIA	Diamond Phototherapy Sheet for Incubator			
HC1MINC002UBIA	Incubator Sheet	75 x 39.5 x 6 cm; 29.5" x 15.5" x 2.3"	Dräger Isolette™ 8000	100 per box
HC2MINC002UBIA	Diamond Phototherapy Sheet for Incubator			
HC1MINC003UBIA	Incubator Sheet	74 x 36 x 6 cm; 29.1" x 14.2" x 2.3"	ATOM V2100-G; ATOM V2200	100 per box
HC2MINC003UBIA	Diamond Phototherapy Sheet for Incubator			
HC1MINC004UBIA	Incubator Sheet	61 x 40 x 6 cm; 24" x 15.7" x 2.3"	NATUS NatalCare™ LX	100 per box
HC2MINC004UBIA	Diamond Phototherapy Sheet for Incubator			
HC1MINC005UBIA	Incubator Sheet	63.5 x 35 x 6 cm; 25" x 13.8" x 2.3"	NATUS NatalCare™ PC-305	100 per box
HC2MINC005UBIA	Diamond Phototherapy Sheet for Incubator			
HC1MBAS001UBIA	Bassinet Sheet	66 x 31.5 x 6 cm; 26" x 12.4" x 2.3"	Most Hospital Bassinets	100 per box
HC1SPM004UBIA	Crib Sheet	85 x 68 x 10 cm; 33.5" x 26.8" x 4"	Most Hospital Cribs	100 per box
HC1SPM005UBIA	Crib Sheet	85 x 59 x 10 cm; 33.5" x 23.2" x 4"	Most Hospital Cribs	100 per box

## Selling Policy

1. Minimum order is \$300.00. For orders under \$300.00, a \$20.00 processing charge will be added to the invoice.
2. All shipments are F.O.B. Buffalo Grove, IL 60089, or point of origin.
3. Payment Terms are net 30 days.
4. Freight will be prepaid on orders of \$5,000.00 or more when shipped to a single location within the continental U.S.A. Shipments less than \$5,000.00 will be shipped freight collect.
5. A fee will be added to the invoice for drop shipments:
  - a. For orders up to \$400.00—\$20.00
  - b. For orders over \$400.00—5% of invoice total

## Returned Goods Policy

It is understood that the responsibility for inventory levels and rotation of inventories rests with the distributor and not with CORPAK MedSystems Corporate office. All returns must have freight prepaid.

1. RMAs for full credit will be issued only for defective merchandise or orders processed in error by CORPAK MedSystems. These returns must be made within 30 days of the invoice date.
2. Any other return requests are subject to the following conditions:
  - a. A 25% restocking fee will be charged for goods returned within 60 days of invoice date. Merchandise must be in salable condition.
  - b. No credit will be issued for out of date goods, modified or discontinued goods. Do not return these goods to CORPAK MedSystems.
  - c. No credit will be issued for goods over 90 days old from the invoice date.
  - d. Return request for dated items are granted at the discretion of CORPAK MedSystems.
3. ENFit™ (ISO Standard 80369-3) Addendum
  - a. No returns on legacy products, other than for defective merchandise or orders processed in error by CORPAK MedSystems within 30 days of the invoice date. Legacy products include all products manufactured without the new ENFit™ Enteral access connector.
  - b. ENFit™ transition connectors are available in the market and provide connection between ENFit™ enabled administration sets and non-ENFit™ feeding tube connections – CORPAK legacy feeding tube products.

## Your Clear Connection to Enteral Safety: Reducing the Risk of Misconnections

CORPAK MedSystems is proud to continue its legacy of leadership in bringing Enteral safety solutions to market. As the innovator for the market leading CORFLO yellow feeding tubes and the CORTRAK Enteral Access System (real-time, non-radiographic visualization technology for feeding tube placement), CORPAK knows how important safe Enteral therapy is to the clinicians and patients we serve every day.

## A New Standard for Enteral Safety

The International Organization for Standardization (ISO) established ISO standard 80369-1 as the framework for all small bore connectors to reduce medical tubing misconnections. Most recently, an international group of clinicians, manufacturers, suppliers, and the U.S. Food and Drug Administration (FDA) have partnered with the ISO and the Association for the Advancement of Medical Instrumentation (AAMI) to create an ISO standard specifically for Enteral connectors (80369-3). The new Enteral connector for the new ISO standard (80369-3) has been given the name “ENFIT”. CORPAK has established a leadership role with industry on this initiative as a founding and charter member of the Global Enteral Device Supplier Association (GEDSA).

## Learning Resource and Training Materials



CORPAK University  
www.corpakU.com


CORPAK, CORPAK logo, CORPAK University, CORPAK University logo, CORTRAK, CORTRAK logo, EAS, CORFLO, CORFLO logo, CORLITE, CORLITE logo, CORGRIP, CORGRIP logo, ANTI-IV, CORLOCK, CORPORT, CORSCOPE, CORFLO-cuBBy, NAVIGATOR, NAVIGATOR logo, MAPCATH, FARRELL, FARRELL logo, and C-19 are either registered trademarks or trademarks of CORPAK MedSystems, Inc. in the United States and/or other countries.

TubeClear® is a trademark of Actuated Medical, Inc., Lifenes® is a trademark of Ubimed, Inc., Skin-Prep™ is a trademark of Smith & Nephew, Inc., Giraffe™ is a trademark of General Electric, Isolette™ and Caleo™ are trademarks of Dräger Medical System, Inc., Natus NatalCare™ is a trademark of Natus Medical, Inc., and ENFit™ is a trademark of GEDSA.





# Your **CLEAR** Connection to Enteral Safety



**CLEAR**  
visualization of  
the fluid path  
and connection

Multiple Patents Pending

[ENFITsafety.com](http://ENFITsafety.com)

1001 Asbury Drive • Buffalo Grove, IL, 60089 • +1 847-403-3400 • Toll Free: +1 800-323-6305 • Fax: +1 847-541-9526

[www.corpak.com](http://www.corpak.com)